

Informatics in Radiology

RADIANCE: An Automated, Enterprise-wide Solution for Archiving and Reporting CT Radiation Dose Estimates¹

TEACHING POINTS

See last page

Tessa S. Cook, MD, PhD • Stefan L. Zimmerman, MD • Scott R. Steingall, RT
Andrew D. A. Maidment, PhD • Woojin Kim, MD • William W. Boonn, MD

There is growing interest in the ability to monitor, track, and report exposure to radiation from medical imaging. Historically, however, dose information has been stored on an image-based dose sheet, an arrangement that precludes widespread indexing. Although scanner manufacturers are beginning to include dose-related parameters in the Digital Imaging and Communications in Medicine (DICOM) headers of imaging studies, there remains a vast repository of retrospective computed tomographic (CT) data with image-based dose sheets. Consequently, it is difficult for imaging centers to monitor their dose estimates or participate in the American College of Radiology (ACR) Dose Index Registry. An automated extraction software pipeline known as Radiation Dose Intelligent Analytics for CT Examinations (RADIANCE) has been designed that quickly and accurately parses CT dose sheets to extract and archive dose-related parameters. Optical character recognition of information in the dose sheet leads to creation of a text file, which along with the DICOM study header is parsed to extract dose-related data. The data are then stored in a relational database that can be queried for dose monitoring and report creation. RADIANCE allows efficient dose analysis of CT examinations and more effective education of technologists, radiologists, and referring physicians regarding patient exposure to radiation at CT. RADIANCE also allows compliance with the ACR's dose reporting guidelines and greater awareness of patient radiation dose, ultimately resulting in improved patient care and treatment.

©RSNA, 2011 • radiographics.rsna.org

Abbreviations: ACR = American College of Radiology, CTDI = CT dose index, DICOM = Digital Imaging and Communications in Medicine, DIR = Dose Index Registry, DLP = dose-length product, IHE = Integrating the Healthcare Enterprise, NRDR = National Radiology Data Registry, OCR = optical character recognition, PACS = picture archiving and communication system, PC = personal computer, RADIANCE = Radiation Dose Intelligent Analytics for CT Examinations, RDSR = radiation dose structured report, REM = Radiation Exposure Monitoring

RadioGraphics 2011; 31:1833–1846 • Published online 10.1148/rg.317115048 • Content Codes: **CT** **IN** **QA**

¹From the Department of Radiology, Hospital of the University of Pennsylvania, 3400 Spruce St, 1 Silverstein, Philadelphia, PA 19104 (T.S.C., S.R.S., A.D.A.M., W.K., W.W.B.); and Department of Radiology, Johns Hopkins University, Baltimore, Md (S.L.Z.). Presented as an education exhibit at the 2010 RSNA Annual Meeting. Received March 11, 2011; revision requested May 26 and received July 11; accepted July 28. T.S.C., S.R.S., A.D.A.M., W.K., and W.W.B. have disclosed various financial relationships (see p 1845); the other author has no financial relationships to disclose. **Address correspondence to** T.S.C. (e-mail: tessa.cook@uphs.upenn.edu).

Funding: The work was supported in part by the National Institutes of Health [grant number T32-EB-004311-06].

Introduction

The use of computed tomography (CT) has increased significantly in the past decade (1). Patients presenting to the emergency department often undergo CT as part of their work-up, and patients with malignancies or other chronic disease may undergo annual or more frequent follow-up imaging (2). In the past 2 decades, the proportion of background radiation in the United States attributed to medical imaging has increased from approximately 15% in 1987 to nearly 50% today (3,4). With the increasing application of CT, ranging from screening to cancer follow-up to emergent evaluation of acute coronary syndrome, patients are being irradiated more often and from an earlier age. The current generation will likely receive a higher lifetime radiation dose than did their predecessors, and the risks associated with this increased dose are difficult to estimate (5,6).

As the lay press and state and federal legislatures demonstrate growing interest in the ability to monitor and report the dose associated with imaging-related radiation (7–9), professional societies have expanded their efforts to educate both patients and physicians on this topic. The ImageGently and ImageWisely campaigns of the Society of Pediatric Radiology and the American College of Radiology (ACR), respectively, are raising awareness in both the medical community and the patient community at large regarding diagnostic irradiation (10,11). The American Association of Physicists in Medicine is working on methods and guidelines for standardization of CT protocols to optimize dose reduction (12). Recent legislation in the state of California requires that radiologists include dose information in their CT reports beginning in July 2012 (13).

Technical developments over the past decade have allowed greater standardization of dose reporting. The Digital Imaging and Communications in Medicine (DICOM) Structured Reporting standard defines dose objects dedicated to storing CT radiation dose information (14,15). The dose objects contain information about the study type, in addition to the relevant CT physics parameters for individual irradiation events (discussed in the next section). The Integrating the Healthcare Enterprise (IHE) initiative has developed a Radiation Exposure Monitoring (REM) profile to assist vendors in increasing the interoperability of their systems with respect to dose monitoring and reporting (16). IHE REM defines a set of “actors” that produce, report, and

transmit dose information, including the acquisition modality, image manager-archive, dose information consumer, dose information reporter, and dose registry. It also defines the interactions between the actors, including DICOM query-retrieve operations and secure transfer of dose information to a registry, so that vendors can develop software and hardware that complies with IHE REM and communicates seamlessly with other compliant systems.

In addition, large-scale initiatives are underway to promote dose monitoring and dose reporting. The ACR Dose Index Registry (DIR), part of the National Radiology Data Registry (NRDR) and often referred to as the NRDR-DIR, is currently collecting population-based data from facilities across the United States (17,18). These data will allow development of reference dose levels and comparison of dose estimates among imaging centers, thereby fostering development of dose optimization measures. The Initiative to Reduce Unnecessary Radiation Exposure from Medical Imaging was recently launched by the U. S. Food and Drug Administration (19). The NIH (National Institutes of Health) is also making efforts to track and report radiation doses for all patients imaged at the Institutes (20).

Despite these various endeavors to initiate more careful monitoring of patient radiation dose in the future, every imaging center has a vast repository of retrospective CT data that does not currently include such information in an easily accessible form. Dose-related parameters for CT have historically been stored as an image-based DICOM dose report, rather than as structured data within the DICOM header. As a result, it is difficult for imaging centers to participate in the ACR DIR. Furthermore, CT scanners currently in use may not have firmware amenable to incorporating radiation dose into image headers. At our institution alone, there are more than 800,000 CT examinations stored in our picture archiving and communication system (PACS) that date back to the early 2000s for which radiation dose data exist but are not easily accessed. In addition, imaging centers without plans to purchase new scanners in the near future will face the challenge of how to most efficiently and accurately record and report radiation dose information.

In this article, in the interests of improving patient care by decreasing exposure to unnecessary medical radiation, improving the benefit-to-risk ratio for patients undergoing diagnostic imaging, and addressing the growing awareness of potential problems caused by radiation exposure at CT,

we describe Radiation Dose Intelligent Analytics for CT Examinations (RADIANCE), an automated extraction software pipeline that we developed for the extraction, archiving, and monitoring of radiation dose information for CT.

CT Parameters That Affect Dose

Although CT dose parameters are regularly reported when a CT scan is performed, either at the CT scanner console or subsequently in the form of an image-based dose sheet, an institutional review board–approved survey of radiologists and radiology trainees conducted at our institution revealed that these data are not routinely reviewed at the time of study interpretation. Despite the considerable variation among vendors with respect to the format of the CT dose sheet used, enough information is typically provided to allow estimation of patient dose for a particular examination. Among physicians, radiologists receive the most training in imaging physics and are thus best equipped to interpret the dose parameters and their effects on image quality as well as patient dose. In this section, we review the relevant dose parameters and important concepts in CT physics.

The kilovolt peak represents the maximum voltage applied across the x-ray tube during imaging and is an important determinant of radiation dose. The kilovolt peak determines the energy of the x-ray beam: The higher the kilovolt peak, the higher the energy of the beam. Altering the kilovolt peak affects image quality: Higher-energy photons will penetrate the patient and contribute to the image, whereas lower-energy photons will be absorbed in the patient before reaching the detector. Because dose is proportional to the kilovolt peak squared, relatively small changes in kilovolt peak can result in large perturbations in dose. Most routine adult examinations use a kilovolt peak of 120; however, imaging at a lower kilovolt peak (100 or even 80 kVp) has been implemented in thinner adults and can greatly reduce dose. Studies performed with low-kilovolt-peak imaging have demonstrated dose reductions of up to 45%–50% for patients undergoing cardiac CT (21,22). Lower-kilovolt-peak imaging is particularly advantageous in vascular imaging because the increased iodine absorption achieved with the lower-energy beam results in higher attenuation of contrast material–enhanced vessels. However, the tradeoff is increased image noise, particularly in large patients, as well as more blooming and streak artifacts from metal or calcification. Artifacts from calcification become problematic in the evaluation of heavily calcified coronary arteries or coronary artery stents.

Milliamperere-seconds, also known as the tube current–time product, determine the number of photons generated by the scanner over a given scanning time. Dose varies linearly with tube current. Fewer milliamperere-seconds will always decrease patient dose at the expense of increased image noise, since fewer photons reach the x-ray detectors. Noise increases by a factor of 1 divided by the square root of total milliamperere-seconds, so that reducing milliamperere-seconds by one-half will increase noise by approximately 40%. Most modern CT scanners use automatic tube current modulation to reduce patient dose (23). Tube current modulation can modify tube current in one of two ways. Along the z-axis of the patient (ie, from head to toe), increasing the number of photons when imaging higher-attenuation regions of the body (eg, shoulders) and decreasing the number of photons when imaging lower-attenuation regions (eg, lungs) can reduce overall patient dose. In addition, angular tube current modulation (in the x-y plane) modifies tube current based on the angle of the x-ray source relative to the rotation arc. Tube current increases when the beam is parallel to the scanner table and must take a long path through the body. It is reduced when the beam is perpendicular to the scanner table and travels a shorter distance to penetrate the patient. These automated techniques can result in dose reductions of 20%–40% (22).

The CT dose index (CTDI), expressed in milligray, indicates the energy output ($1 \text{ Gy} = 1 \text{ J/kg}$) of the CT scanner in different situations. CTDI_{100} is the most straightforward and refers to the energy imparted to a 100-mm ion chamber placed in the center of a standard cylindric acrylic phantom during a CT examination. CTDI_w refers to a weighted average of energy output to the entire phantom, taking into account the incident energy differences between the periphery (higher energy) and center (lower energy) of the phantom. Finally, CTDI_{vol} represents a further refinement of CTDI_w that accounts for pitch. CTDI_{vol} is influenced by kilovolt peak and milliamperere-seconds in addition to pitch. **It is important to remember that CTDI_{vol} is not a measure of patient dose, but rather represents the energy output of the scanner (24).**

The dose-length product (DLP), expressed in milligray-centimeters, represents the dose to the cylindric CT phantom and can be used to estimate patient dose at CT. The DLP is the scan length (in centimeters) multiplied by CTDI_{vol} . It

can be converted to an estimated patient effective dose (in millisieverts) by multiplying by unitless anatomy-specific conversion factors (also known as *k* factors) (25). These conversion factors are derived from tissue-weighting factors that are maintained by the International Commission for Radiation Protection (26) and evolve as the Commission updates the factors. Studies that irradiate the thyroid gland, breasts, or gonads typically have higher *k* factors than do those that irradiate the brain or extremities. The change in units for effective dose from milligray to millisieverts reflects the application of the tissue-weighting factors. With use of the DLP, it is possible to derive an estimated effective dose ($DLP \times k$, or $CTDI_{vol} \times \text{scan length} \times k$) for a patient's CT examination. Because DLP is dependent on scan length, reducing the scan length will lower the patient dose but may also introduce the risk of missing relevant disease by excluding it from the imaged anatomy.

Additional nonphysics parameters, such as patient age and patient size, can also affect dose (27,28). Younger patients are at a higher risk of experiencing the adverse effects of radiation exposure because of more rapid turnover of cellular DNA and the higher potential for radiation-induced mutations. Similarly, women who undergo breast imaging are more susceptible. In addition, for a given CT protocol, a thinner patient will receive a higher organ dose than will a larger patient owing to a smaller amount of attenuating soft tissue. Hence, CT protocols should always be tailored to patient size so as to (a) optimize diagnostic image quality and dose savings, and (b) maximize the benefit-risk ratio of exposing patients to ionizing radiation.

RADIANCE Algorithm

Figure 1 illustrates the RADIANCE pipeline. The input to the pipeline is the DICOM CT dose sheet. Scanner manufacturers typically store these documents as a separate series within the study, thereby simplifying their retrieval from a PACS. Some vendors use a unique series number; for example, series 501 for Siemens Medical Systems (Erlangen, Germany), series 999 for GE Medical Systems (Milwaukee, Wis), and (more recently) series 9000 for Toshiba Medical Systems (Tustin, Calif). Neither dose sheets generated by Philips Medical Systems (Best, the Netherlands) nor Toshiba dose sheets generated before 2009 are stored in series that are consistently numbered.




Figure 1. Chart illustrates the RADIANCE automated dose extraction pipeline, which combines data from the dose sheet, examination header, and radiology information system (RIS) to allow analytics and quality assurance. *OCR* = optical character recognition, *XML* = extensible markup language.

The first step in the pipeline is optical character recognition (OCR) analysis of the dose sheet. This step converts the pixel-based information into ASCII text (29), which can be parsed to extract the relevant CT dose parameters. The OCR tool also translates the dose sheet's DICOM header into extensible markup language and produces a JPEG (Joint Photographic Experts Group) copy of the original DICOM dose sheet for validation. These copies are used to validate the performance of the software (30). For some vendors (eg, Philips and Toshiba), a binary threshold must be applied to the dose sheet before OCR is applied, since the dose data are stored in the actual pixels of the image instead of in an overlay, as with Siemens.

The relational database stores an aggregate of study data by parsing both the dose sheet and the study header. Before the OCR version of the dose sheet is parsed, useful study parameters are parsed from the extensible markup language study header to tailor the processing steps that follow. The information extracted from the study header includes the study date, examination code and description, scanner manufacturer and model, institution name, patient identification number, and patient date of birth. This information can be used for quality control and correct insertion of the dose-related values into the database.

Each dose sheet is processed in a vendor-specific fashion (discussed later). The extracted dose-related parameters for each series of a CT examination include kilovoltage, x-ray tube milliamperes, reference milliamperes, $CTDI_{vol}$, and DLP. Each vendor's dose sheet reports a subset

04-Sep-2010 11:05

Ward: OTU
Physician: VASCULAR
Operator: ED

Total mAs 2449 Total DLP 151 mGy*cm

	Scan	kV	mAs / ref.	CTDIvol mGy	DLP mGy*cm	TI s	cSL mm
Patient Position F-SP							
Topogram	1	120	50 mA			5.3	0.6
DS_CaScSeq	2D	120	69 / 76	4.00	50	0.2	3.0
Last scan no.	8						
PreMonitoring	9	100	40	1.07	1	0.33	10.0
I.V. Bolus							
Monitoring	10	100	40	5.33	5	0.33	10.0
Last scan no.	14						
DS_CorCTAAdapt	15D	100	212 / 324	7.97	95	0.2	0.6
Last scan no.	21						

Figure 2. Dose sheet from a coronary CT examination performed on a Siemens scanner. Siemens dose sheets contain no protected health information.

of these parameters. During software development, consistently repeated errors in the OCR process were iteratively identified and corrections were implemented in the code on a vendor-specific basis. Once the provided dose parameters are parsed, validated, and stored in the database, the estimated whole-body effective dose is calculated by multiplying the total study DLP by the appropriate anatomy-specific k factor. RADIANCE uses the k factors derived from International Commission on Radiologic Protection publication 60 (26).

A radiation dose structured report (RDSR) can be produced for every study with an image-based dose sheet. This is achieved using a Java utility built on the PixelMed Java DICOM Toolkit (31). **RADIANCE populates the structured report template with the data from its database and generates an RDSR that can be automatically transmitted to the ACR DIR.** This enables imaging facilities to automate transfer of their dose data to the DIR, eliminating the need for manual data entry or manual transmission of dose reports. **The contents of scanner-generated RDSRs can also be imported into the database, enabling users to centralize their dose monitoring in cases in which there is a mixture of CT scanners, only some of which can produce structured reports.**

Vendor-specific Implementation Details

The parameters extracted from the study header guide the subsequent parsing of the OCR output dose sheet. **RADIANCE is able to process dose sheets generated by the four major CT vendors—Siemens, GE, Toshiba, and Philips—as well as from the CereTom portable head CT scanner by NeuroLogica (Danvers, Mass).** Each vendor's dose sheets have certain unique features that re-

quire customized processing. These specifications are described in the sections that follow.

Siemens

Of the dose sheets generated by the four major vendors, the Siemens dose sheet (Fig 2) is the only one that does not include any protected health information. The dose sheet is found in series 501 of the CT examination (“Patient Protocol”) and typically consists of only one page. In rare instances, multipage dose sheets are generated. The first few lines of the Siemens dose sheet typically include the study date and time as well as information pertaining to the ordering physician, ordering department, and performing technologist. The latter field (“Operator”) usually contains the initials of the technologist who initiated the study on the scanner, but it can store full names as well. The data pertaining to the ordering physician and department are variable; for example, the “Physician” field may include the name of either an individual or a department. Similarly, the “Ward” field usually contains an abbreviation rather than the full name of a department. RADIANCE ignores these values, instead obtaining this information from the radiology information system, where it is more accurate.

The total milliamperere-seconds and total DLP are often included at the top of the dose sheet, although their position can vary with the scanner model. The accession number is typically not provided on the dose sheet.

The remainder of the dose sheet is in a tabular format and includes the scan number, kilovoltage, tube current milliamperes, reference milliamperes, CTDI_{vol}, DLP, rotation time, and

Teaching Point

Teaching Point

Figure 3. Dose sheet from a chest CT angiographic examination performed on a GE scanner. Identifiable patient information that is normally included on the dose sheet has been removed for compliance with the Health Insurance Portability and Accountability Act.

Patient Name:			Exam no:		
Accession Number:			Sep 22 2010		
Patient ID:			LightSpeed VCT		
Exam Description: CT Chest Angiography					
Dose Report					
Series	Type	Scan Range (mm)	CTDIvol (mGy)	DLP (mGy-cm)	Phantom cm
1	Scout	-	-	-	-
200	Axial	199.250-199.250	11.81	5.90	Body 32
2	Helical	511.750-1244.500	23.43	708.99	Body 32
Total Exam DLP:				714.89	

collimated section length for each of the scout images and series within the study. If both tube current milliamperes and reference milliamperes are specified, they are usually delimited by a forward slash. If reference milliamperes are not provided (ie, dose modulation is not turned on), the forward slash is omitted and only one value is provided.

To parse the Siemens dose sheet, RADIANCE first parses the header to extract the operator's identity, total milliamperere-seconds, and (when specified) total DLP. RADIANCE looks for the word "Patient" or "Scan" to locate the start of the dose parameters. Each line of dose information is first examined to see if it matches one of a set of prescribed patterns. These patterns include scout or topographic image entries, which report kilovoltage and (rarely) milliamperes but not CTDI_{vol} or DLP. In addition, special dual-energy sequences have a unique format in which there is a series name for one of the kilovolt settings but not for the other. In addition, the presence of a value for reference milliamperes is specifically sought. Any series that consists only of a series name and number but no dose parameters is omitted. Series that don't fit any of the predefined patterns are parsed according to a standard template of expected parameters. Some versions of the Siemens dose sheet now specify the phantom type used for the study by providing a footnote of "(a)" or "(b)" or, more recently, "S" or "L" with each CTDI_{vol} and defining the footnote at the bottom of the dose sheet. If this information is specified, it is also stored in the database.

The volume and rate of intravenous contrast material injection are sometimes reported on the dose sheet, either as one of the series entries with the dose information or separately at the bottom of the dose sheet. At present, this information is not stored in RADIANCE.

Patient Name (Country) :

Patient Name (Multi-byte) :

```

ID :                               Study ID :
Birth Date : 1961. .              Age : 48Y
Sex : F   Weight(kg) :             Height(cm) :
Patient Comments : PANCREATITIS
Study Date : 2010.04.20           Body Part : ABDOMEN/PELVIS
Requesting Department :
Referring Physician :
Reporting Physician :
Operator Name :
Total Image Number : 1384
Image Comment :

<< Dose Information >>
Total mAs : 7842                    Total Scan time : 45.38
CTDIvol (Head) : - (Body) : 1285.10
DLP (mGycm) (Head) : - (Body) : 1001.50


<< Contrast/Enhance Information >>
Contrast Enhance : CE

```

Figure 4. Dose sheet from an abdominopelvic CT examination performed on a Toshiba scanner. Protected health information that is normally included on the dose sheet has been removed for Health Insurance Portability and Accountability Act compliance.

GE Medical Systems

The GE dose sheet (Fig 3) is found in series 999 ("Dose Report"). The header section of the dose sheet often includes the patient's name and identification number, as well as the study accession number, examination description, scanner model name, and study date. RADIANCE does not parse any of these data from the dose sheet because all of this information is included in the DICOM header, where it is immune to potential errors in the OCR process. The "Operator" field specified on the Siemens dose sheet is not included by GE and defaults to a value of "ZZZ" to indicate that these data were not specified on the



#	Description	Mode	CTDI [mGy]	DLP [mGy*cm]
2	SURVIEW	Surview	0.0	0.00
3	CHEST, ABD, PEL	Helical	11.6	839.00
3	CHEST, ABD, PEL	Helical	17.6	220.88
6	FEMUR	Helical	16.5	1149.70

Figure 5. Dose sheet from a chest and abdominopelvic CT examination performed on a Philips scanner. Aside from the study identification number, Philips dose sheets do not contain any protected health information.

dose sheet. In addition, total milliampereseconds are not specified and default to an unspecified value of -1 .

RADIANCE searches for the term “mGy” to locate the start of the dose parameters. The remainder of the GE dose sheet is tabular and reports a number, type, scan length, $CTDI_{vol}$, DLP, and phantom type for every series in the study. Of these data, RADIANCE stores the series type (ie, series name), $CTDI_{vol}$, DLP, and phantom type. At present, the scan length is not included in the database, since its value is inherently reflected in the series DLP. Neither tube current nor reference tube current is specified on the dose sheet, and they default to -1 in the database. Unlike the Siemens dose sheet, on which not all series entries follow the same tabular format, the GE dose sheet has no such discrepancy.

The last line of the GE dose sheet reports the total DLP for the examination. This value is also parsed and stored. It can be compared with the sum of the individual series DLPs for internal validation of the OCR process as well as the parsing algorithm.

In some cases, GE studies produce multipage dose sheets. In such cases, the number at the bottom of each page of the dose sheet will specify the total number of pages in the dose sheet. All of the individual pages will be found within series 999 of the study.

Toshiba

Historically, the location of the Toshiba dose sheet (Fig 4) in a CT examination has varied. For studies from approximately 2009 onward, they can be found in series 9000 of the examination, but they may not have a consistent series name. For earlier examinations, both the series name and the series number can vary. This poses a challenge in retrieving these dose sheets from the PACS.

Parsing of the dose sheets is relatively straightforward. The Toshiba dose sheet headers are quite detailed, but as with the GE dose sheets, most of the data are also available in the DICOM header, which is not susceptible to OCR extrac-

tion errors. However, the operator’s name or initials are parsed from the dose sheet and stored in the database.

The total milliampereseconds, average $CTDI_{vol}$, and total DLP for the study, in addition to the type of phantom on which these estimates were based (ie, head or body), are parsed from the first page of the dose sheet. In some cases, the sum of the $CTDI_{vol}$ values for individual series is provided; although these data cannot be used to correctly calculate DLP for each body part scanned, they are extracted and stored. Users with these scanner models should verify how the $CTDI_{vol}$ is reported to ensure correct use of the data for analysis. As with other vendors, information related to the use of intravenous contrast material is not parsed from the dose sheet. Reference tube current is not provided. Newer Toshiba dose sheets have a multipage format in which the first page retains the original summary format and subsequent pages provide series-specific values.

Philips

Dose sheets from Philips scanners (Fig 5) are found in a series named “Dose Info” but are not always consistently numbered from one study to the next. Aside from the study identification number, the header of the Philips dose sheet contains no identifiable patient information; the information under the header is in a very consistent tabular format. One unique feature of the Philips dose sheet headers is the “Estimated Dose Savings” field, which is, however, currently not archived within RADIANCE.

The total DLP is parsed from the Philips dose sheet header, but the additional information, which includes the study identification number, accession number, and study date and time, is obtained from the DICOM study header. RADIANCE looks for the term “CTDI” to locate the start of the individual series’ dose parameters. Only the series name and number, acquisition mode (eg, scout, axial, or helical), $CTDI_{vol}$, and DLP for each series are reported on the dose

sheet; no phantom information is provided. For each series in the study, the software archives the series name, CTDI_{vol}, and DLP.

Another unique feature of Philips scanners is that they have historically provided dose parameters in the DICOM header, although not in the same format as the new RDSR. Nevertheless, these data are also parsed and are used to validate the values extracted from the dose sheet by means of OCR.

NeuroLogica

The dose sheets from a CereTom portable head CT scanner (NeuroLogica) (Fig 6) are not consistently numbered. They are typically found in either series 3 or series 4 within the PACS, in a series named "Dose Report." RADIANCE looks for the term "Range" to locate the start of the series-related dose parameters. The dose sheet reports the number, type (eg, axial), scanning range, CTDI_{vol}, and DLP for each series. In addition, the total DLP is reported on a separate line. RADIANCE archives the series type, CTDI_{vol}, series DLP, and total DLP. Kilovoltage, tube current milliamperes, reference milliamperes, and phantom type are not reported on CereTom dose sheets; instead, the numeric data are stored with default unspecified values of -1 and the phantom type is stored as "Unknown."

RADIANCE Technical Details

Implementation

RADIANCE requires installation of a World Wide Web server and a MySQL database server to run; development was performed using XAMPP, an open-source software that provides this functionality (32). The software is implemented using the PHP scripting language, which is platform independent. Windows binaries (Microsoft, Redmond, Wash) for the OCR software are used, and the automated pipeline script is a Windows batch file, which in future versions will also be implemented for Linux (Linus Torvalds, Santa Clara, Calif) and Macintosh (Apple Computer, Cupertino, Calif). The pipeline batch file can be scheduled as a recurrent Windows task (using the "Scheduled Tasks" feature in Windows) to run as frequently as desired by the imaging facility. Independent of the OCR step, RADIANCE runs successfully on both Windows and Macintosh platforms but has yet to be tested on Linux. However, there are no inherent platform dependencies within the PHP code.

Dose Report					
Series	Type	Scan Range (mm)	CTDIvol (mGy)	DLP (mGy-cm)	
1	Axial	-382.04 - -192.04	47.15	895.89	
				Total:	895.89

Figure 6. Dose sheet from a portable head CT examination performed on a NeuroLogica CereTom scanner. CereTom dose sheets do not contain any protected health information.

RADIANCE has a small footprint and can run on any standard, up-to-date Windows personal computer (PC). This PC must be configured as a DICOM destination, that is, assigned an application entity title with an open port to receive DICOM dose sheets directly from scanners or from a PACS. The scanner or PACS must be able to recognize and communicate with this application entity title; a onetime manual authentication may need to be configured. Dose sheets can be queried by the PC or sent from the PACS to the RADIANCE inbox, a folder within the installation. This inbox is searched for new DICOM dose sheets, which are then processed by means of the OCR software. The OCR output and the DICOM study header are then automatically parsed, and extracted data are stored in the database. This implementation allows automatic, behind-the-scenes dose monitoring that does not interfere with the technologist's or radiologist's work flow.

To perform DICOM query-retrieve operations from a PC running RADIANCE, additional software is needed to either (a) run the PC as a DICOM listener (using the *dcm4che2* toolkit [33]), or (b) query and retrieve dose sheets directly from the PACS (using the *dcm4che2* toolkit [34]). A DICOM listener is easily configured using the *dcm4che2* utility within *dcm4che2*, which allows assignment of an application entity title and listens for data (in this case, DICOM dose sheets) on an open port. *dcm4che2* provides utilities for standard DICOM C-FIND and C-MOVE operations (*findscu* and *movescu*, respectively) to allow both study- and series-level DICOM queries, identification of specific series numbers, and retrieval of dose sheets, which becomes challenging when dose sheets are not consistently named or numbered.

During the IHE Dose Demonstration at the 2010 Annual Meeting of the Radiological Society of North America, RADIANCE was shown to be compliant with the IHE REM profile. Within the profile, the software assumed the roles of two actors: dose information consumer and dose information reporter (described in greater detail later).

Teaching
Point




Figure 7. RADIANCE database schema. Each box represents a table within the database. The primary key for each table is shown in bold. The *masterflatfile* table, a combination of all the data points in the *masterscanner*, *masterprovider*, *masterpatient*, and *mastertotals* tables, is not shown.

Database Schema

The RADIANCE core database schema is shown in Figure 7. The study accession number serves as the primary key for most tables in the database. Because the study accession number appears multiple times in the *masterdose* table (once in association with each series of a study), an automatically incremented primary key is used instead (Fig 7).

The tables that constitute the primary functionality of RADIANCE are named with the prefix “master.” The *masterdose* table stores the dose-related parameters for the series of a CT examination, each of which is identified by its accession number. The *masterpatient* table records one entry for each accession number, rather than for each patient. It includes the patient’s height and weight, which can change from one examination to the next. As a result, a unique entry is made for each study the patient undergoes. The *master-scanner* table includes information about the specific scanner on which the study was performed, and notes both the manufacturer and the model name as well as the station name, or unique identifier, for a scanner at a particular institution. The *masterprovider* table is populated predominantly by querying the departmental radiology informa-

tion system for information on the personnel who requested, performed, and interpreted the study. The operator’s name or initials, which are sometimes specified in the DICOM study header, are also stored here.

The *mastertotals* table includes the total milli-ampere-seconds (if provided) and the total DLP, both as reported (if available) and as calculated (ie, summed from the series in the study). The anatomy-specific conversion factor used to estimate whole-body effective dose, as well as this dose estimate, are also stored in the table. The body parts imaged in a study are matched up with the corresponding conversion factor with the use of two lookup tables: *examcodeslookup* and *publisheddoselookup*.

A *masterflatfile* table is constructed by combining the contents of the *masterscanner*, *master-provider*, *masterpatient*, and *mastertotals* tables, since all of these tables share the accession number of the study as a field. The *masterflatfile* table is suitable for exporting data from the database for subsequent analysis and is also used by the dashboard (see “Additional Features of RADIANCE”).

Figure 8. Validation tool used to assess the accuracy of RADIANCE and identify incorrect parsing of dose sheet parameters. Top right-hand panel shows the dose-related parameters extracted for each series within the study, as well as the total DLP (when provided). Bottom right-hand panel allows the user to identify common errors, such as a mismatch between reported and total DLP, missing series data, and incorrect extraction of the kilovoltage and CTDI. The “Comments” field was provided for reporting uncommon errors. Left-hand panel displays the original dose sheet for comparison.

Validator

1799517

Original Image

04-Aug-2006 14:16

Ward: Physician: POP BODY CBBP Operator:

Total mAs 16066 Total DLP 2528

Scan	kV	mAs / ref	CTDIvol	DLP	TI	cBL
1	120				5.3	0.8
2	120	250	16.85	800	0.5	1.2
3	120	250	19.15	864	0.5	0.8
4	120	250	19.15	864	0.5	0.8

Patient Position F-SP Topogram PRE Contrast NephroPhase EXCRETION

From Database, after OCR

Total mAs	Reported Total DLP	Calculated Total DLP
16066	2528	2528

Series	kV	mAs	RefmAs	CTDI	DLP
Topogram1	120	0	0	0	0
PRE	120	250	0	16.85	800
NephroPhase	120	250	0	19.15	864
EXCRETION	120	250	0	19.15	864

Yes No

Total DLP Correct? ☒ ☐

Calculated DLP Correct? ☒ ☐

Missing Series? ☐ ☒

kV Correct? ☒ ☐

CTDI Correct? ☒ ☐

Comments:

Submit Validator: Record 11 of 100

Additional tables help organize the data that are used to create the RDSRs that are transmitted to the ACR DIR.

RADIANCE Validation

Unique repeatable errors were found during OCR processing of each of the different vendors' dose sheets. For example, the number 7 was always replaced by the letter J when Siemens dose sheets were processed. Similar alphanumeric substitutions were identified for other vendors' dose sheets. All identified errors were consistent and repeatable for a particular vendor; consequently, correction of the parameters extracted from the OCR output was vendor specific. No manual intervention on the part of the user is required to correct these errors; the software expects and accounts for these substitutions.

The accuracy of RADIANCE in extracting dose-related parameters from CT dose sheets was validated by checking the extracted data against the original image-based dose sheet. A validation tool was designed for this purpose (Fig 8).

A total of 2108 studies randomly chosen from our institutional PACS were used in the validation (30). The dose sheets for these studies were processed, after which the extracted total DLP was compared with the total DLP reported on the dose sheets. An accuracy of 100% was achieved during this first phase of validation.

A second phase of validation was performed in which the extraction of additional dose parameters was evaluated. During this phase, errors were detected in the extraction of total milliampereseconds (when reported), as were errors involving parameters for dual-source examinations obtained on Siemens scanners. The OCR parsing algorithm was modified to correctly extract these parameters, resulting in an accuracy of more than 99% over all parameters. Errors most commonly occurred as a result of unexpected shifts in the reading frame used to extract the parameters.

Additional Features of RADIANCE

Compliance with the IHE REM Profile and the ACR's NRDR-DIR

The IHE REM profile is an implementation guide for both vendors and consumers that provides a set of standards for radiation dose monitoring. **RADIANCE is compliant with the IHE REM profile as both a dose information consumer and a dose information reporter.** As a dose information consumer, RADIANCE can be configured to query the PACS for CT dose sheets using standard DICOM query-retrieve operations. It is also able to read RDSRs produced by newer scanner models and import these data into the database. Because the database contains protected health information, it is password protected and is intended to reside

Teaching Point




Figure 9. Screen shot from a patient summary within the RADIANCE dashboard, which provides individual dose estimates as well as a cumulative lifetime patient dose estimate for all imaging undergone at a given facility. This representation of the data gives users a sense of how many and what types of CT examinations a patient has undergone, as well as an indication of who ordered these studies.

behind an imaging center's firewall, rather than on a public network.

As a dose information reporter, the analytics dashboard built on the database schema allows users to scrutinize their dose data more carefully and to analyze dose estimates by departmental section, individual scanner, involved personnel, or individual patient. Outlier identification is also possible, that is, the detection of studies whose dose estimates exceed a prescribed threshold. With the exception of the patient dashboard, which requires input of the patient's medical record number, all data presented by the dashboard are completely deidentified in keeping with the REM profile. One of the actions of a dose information reporter in the REM profile is to transmit dose information to a registry. Hence, RADIANCE is able to generate deidentified RDSR representations of legacy CT dose sheets, thereby enabling users to participate in dose registries. This capacity was tested in the IHE Dose Demonstration at the 2010 Annual Meeting of the Radiological Society of North America, during which RADIANCE fulfilled the roles of both dose information consumer and dose information reporter. With use of secure free-text protocol, RDSR representations of legacy CT dose sheets were produced and transmitted to the ACR DIR to demonstrate communication between a dose information reporter and a dose registry. Alternatively, these representations can be routed

through the ACR's TRIAD Server, installed locally at each facility, and sent to the NRDR-DIR. RDSR generation is implemented using the open-source PixelMed Java DICOM Toolkit (31).

Dashboard Analytics

To facilitate review of the extracted data, we have designed a customizable dashboard built on the RADIANCE database schema. The dashboard is implemented using Highcharts JS (35), an open-source JavaScript charting library (Netscape [<http://en.wikipedia.org/wiki/JavaScript>]). The dashboard provides a set of standard overview screens that summarize dose estimates by departmental section, scanner, personnel, patient, and examination type. In addition, it allows users to identify those examinations that exceed a prescribed threshold for estimated whole-body dose (in millisieverts). Moreover, the dashboard is easily customizable and enables users to add screens to those already provided that organize the data according to the needs at their particular facility.

A screen shot from a patient summary is shown in Figure 9. The patient dashboard indicates the dose estimates for each individual CT examination at a given facility, as well as a cumulative lifetime dose estimate for the patient in question. Despite the limitations of adding dose estimates for separate CT scans, this representation gives




Figure 10. Screen shot from a personnel summary within the RADIANCE dashboard, which displays dose estimates in descending order for referring physicians, performing technologists, and reporting radiologists. The names were randomly generated from an online database to anonymize our institutional radiologists.

the physician a sense, not only of how many examinations a patient has undergone, but also of what types of studies were performed and by whom they were ordered. A screen shot from a personnel overview is shown in Figure 10. The personnel overview allows users to examine dose estimates by referring provider, performing technologist, and reporting radiologist, to look for trends in ordering and reporting as well as compliance with protocol implementations. A screen shot from a scanner overview is shown in Figure 11. This overview is useful for identifying (a) potential protocol differences between scanners for the same study type, or possibly (b) a lack of adherence to prescribed protocols for a particular study.

Open-Source Release

RADIANCE debuted as a free, open-source software package in December 2010. The customizable dashboard was released in January 2011, also as a freely available, open-source package. Both are available online at <http://www.radiance-dose.com>. A Google group exists for the community of users and can be accessed at <http://groups.google.com/group/radiance-users>.

Discussion and Future Work

At present, RADIANCE has some limitations in terms of functionality. RADIANCE version 1.0 is not fully compatible with multipage dose sheets, although this functionality is expected with the next major update. In addition, RADIANCE is limited to the archival of dose-related parameters for CT only. Future work will focus on including other modalities, such as fluoroscopy and radiography. Although the software has been most extensively tested on the Windows operating system, there are no inherent operating system dependencies within the RADIANCE PHP code. To date, auxiliary software such as GOCR has been tested on Windows only. Future efforts will focus on providing Linux compatibility for the OCR process and pipeline scripting.

One of the challenges in correctly estimating radiation dose is being able to correctly identify the body parts that were imaged. RADIANCE expects imaging of separate body regions to be reported on separate dose sheets for correct dose estimation. When a dose sheet includes series from more than one body part (eg, chest, abdomen, and pelvis), the software will correctly estimate the dose but may incorrectly attribute the entire dose estimate to the chest because only one accession number is associated with the study.




Figure 11. Screen shot from a scanner overview within the RADIANCE dashboard, which shows average and maximum dose estimates for a particular study type across multiple scanners at a given facility. This overview is helpful for detecting differences in dose between scanners, which may be due to protocol differences or lack of adherence to prescribed protocols.

We have developed an algorithm that addresses this problem by directly querying the PACS database at our institution and looking for scenarios in which multiple body parts are grouped on a single dose sheet (36). Without additional information from either the facility's PACS or radiology information system, however, the algorithm may not be able to correctly sort these doses by the appropriate body part.

Summary

With a growing need to be able to monitor, archive, and report radiation dose estimates and dose-related parameters from CT examinations, RADIANCE provides a free, open-source solution that can be implemented at imaging facilities worldwide with a small footprint on a standard Windows PC. Its users are able not only to monitor CT dose estimates, optimize protocols, and detect outliers, but also to participate in the ACR DIR, even if their scanners do not produce RDSRs. In addition, at facilities with a mix of older and newer scanners, the software will import RDSRs generated by scanners with updated firmware, enabling facilities to use RADIANCE for internal dose monitoring and reporting regardless of the types of scanners available. RADIANCE is also compliant with the IHE REM profile, fulfilling the roles of both dose information consumer and dose information reporter within the profile. This allows imaging centers worldwide to monitor dose estimates administered at CT, optimize protocols on the basis of a review of the dose data collected, and, ultimately, maximize the benefit-to-risk ratio for patients by reducing unnecessary patient exposure to medical radiation while still performing diagnostic imaging examinations.

Acknowledgments:—The authors especially thank Thomas Chiang, MBA, and Daniel Morton, MSEE, of the Medical Imaging Group in the Department of Radiology at the Hospital of the University of Pennsylvania, who were instrumental in the successful deployment of RADIANCE at Penn, assisting with retrospective retrieval of the CT dose sheets from the PACS and the overall setup of the live application. We also thank Jonathan Haski, BS, and Vishnu Ganesan, BS, for their help in implementing the original version of the validator tool.

Disclosures of Potential Conflicts of Interest.—

A.D.A.M.: *Related financial activities:* none. *Other financial activities:* on the board of Real-Time Tomography and holds patents for anthropomorphic phantoms and octree encoding. **S.R.S.:** *Related financial activities:* none. *Other financial activities:* consultant for the Radiological Society of North America. **T.S.C.:** *Related financial activities:* received research support from the National Institutes of Health. *Other financial activities:* none. **W.K.:** *Related financial activities:* none. *Other financial activities:* on the boards of Montage Healthcare Solutions and iVirtuoso; consultant for Amirsys, Click Diagnostics, and Infiniti Medical; received payment from Terumo Medical. **W.W.B.:** *Related financial activities:* none. *Other financial activities:* on the boards of the Society of Imaging Informatics in Medicine, Montage Healthcare Solutions, and Merge Healthcare; holds stock in Montage Healthcare Solutions and Merge Healthcare.

References

1. Maitino AJ, Levin DC, Parker L, Rao VM, Sunshine JH. Nationwide trends in rates of utilization of noninvasive diagnostic imaging among the Medicare population between 1993 and 1999. *Radiology* 2003;227(1):113–117.
2. Levin DC, Rao VM, Parker L, Frangos AJ, Sunshine JH. Recent trends in utilization rates of abdominal imaging: the relative roles of radiologists and nonradiologist physicians. *J Am Coll Radiol* 2008; 5(6):744–747.

3. Sinclair WK, Adelstein SJ, Carter MW, Harley JH, Moeller DW. Ionizing radiation exposure of the population of the United States: technical report 93. Bethesda, Md: National Council on Radiation Protection and Measurements, 1987.
4. Kase KR. Ionizing radiation exposure of the population of the United States: technical report 160. Bethesda, Md: National Council on Radiation Protection and Measurements, 2009.
5. Brenner DJ, Hall EJ. Computed tomography: an increasing source of radiation exposure. *N Engl J Med* 2007;357(22):2277–2284.
6. Smith-Bindman R, Lipson J, Marcus R, et al. Radiation dose associated with common computed tomography examinations and the associated lifetime attributable risk of cancer. *Arch Intern Med* 2009;169(22):2078–2086.
7. Congressional Subcommittee on Health. Medical radiation: an overview of the issues. Available at: http://energycommerce.house.gov/index.php?option=com_content&view=article&id=1910:medical-radiation-an-overview-of-the-issues&catid=132:subcommittee-on-health&Itemid=72. Accessed March 1, 2010.
8. Bogdanich W, Ruiz RR. FDA to increase oversight of medical radiation. *New York Times*, February 9, 2010.
9. Landro L. Radiation risks prompt push to curb CT scans. *Wall Street Journal*, March 2, 2010.
10. Goske MJ, Applegate KE, Boylan J, et al. The 'Image Gently' campaign: increasing CT radiation dose awareness through a national education and awareness program. *Pediatr Radiol* 2008;38(3): 265–269.
11. Brink JA, Amis ES Jr. Image Wisely: a campaign to increase awareness about adult radiation protection. *Radiology* 2010;257(3):601–602.
12. American Association of Physicists in Medicine. CT protocols. Available at: <http://aapm.org/pubs/CTProtocols/>. Accessed July 1, 2011.
13. Padilla A, Alquist EK. California S.B. no. 1237. Radiation control: health facilities and clinics—records. Available at: http://www.leginfo.ca.gov/pub/09-10/bill/sen/sb_1201-1250/sb_1237_bill_20100929_chaptered.html. Accessed February 28, 2011.
14. DICOM Standards Committee. DICOM standard supplement 127: CT radiation dose reporting. Rosslyn, Va: DICOM Standards Committee, 2007.
15. DICOM Standards Committee. DICOM standard part 16: content mapping resource. Rosslyn, Va: DICOM Standards Committee, 2008.
16. Radiation Exposure Monitoring. Available at: http://wiki.ihe.net/index.php?title=Radiation_Exposure_Management. Accessed March 15, 2010.
17. Amis ES Jr, Butler PF, Applegate KE, et al. American College of Radiology white paper on radiation dose in medicine. *J Am Coll Radiol* 2007;4(5): 272–284.
18. American College of Radiology. National Radiology Data Registry. Available at: http://www.acr.org/SecondaryMainMenuCategories/quality_safety/NRDR.aspx. Accessed March 15, 2010.
19. U.S. Food and Drug Administration. White paper: initiative to reduce unnecessary radiation exposure from medical imaging. Available at: http://www.fda.gov/Radiation-EmittingProducts/RadiationSafety/RadiationDoseReduction/ucm199994.htm#_Toc253092884. Accessed March 15, 2010.
20. Neumann RD, Bluemke DA. Tracking radiation exposure from diagnostic imaging devices at the NIH. *J Am Coll Radiol* 2010;7(2):87–89.
21. Hausleiter J, Meyer T, Hermann F, et al. Estimated radiation dose associated with cardiac CT angiography. *JAMA* 2009;301(5):500–507.
22. Luaces M, Akers S, Litt H. Low kVp imaging for dose reduction in dual-source cardiac CT. *Int J Cardiovasc Imaging* 2009;25(suppl 2):165–175.
23. McCollough CH, Bruesewitz MR, Kofler JM Jr. CT dose reduction and dose management tools: overview of available options. *RadioGraphics* 2006;26(2):503–512.
24. McCollough CH, Leng S, Yu L, Cody DD, Boone JM, McNitt-Gray MF. CT dose index and patient dose: they are not the same thing. *Radiology* 2011;259(2):311–316.
25. Christner JA, Kofler JM, McCollough CH. Estimating effective dose for CT using dose-length product compared with using organ doses: consequences of adopting International Commission on Radiological Protection publication 103 or dual-energy scanning. *AJR Am J Roentgenol* 2010;194(4):881–889.
26. 1990 Recommendations of the International Commission on Radiological Protection. *Ann ICRP* 1991;21(1-3):1–201.
27. Huda W, Vance A. Patient radiation doses from adult and pediatric CT. *AJR Am J Roentgenol* 2007;188(2):540–546.
28. Kalra MK, Maher MM, Toth TL, et al. Strategies for CT radiation dose optimization. *Radiology* 2004;230(3):619–628.
29. GOCR. Available at: <http://jocr.sourceforge.net/>. Accessed March 1, 2010.
30. Cook TS, Zimmerman SL, Kim W, Boonn WW. Automated extraction and reporting of dose information from computed tomography examinations. College Park, Md: American Association of Physicists in Medicine, 2010.
31. Clunie DA. PixelMed Java DICOM Toolkit. Available at: <http://www.dclunie.com/pixelmed/software/>. Accessed March 1, 2011.
32. XAMPP. Available at: <http://www.apachefriends.org/en/xampp.html>. Accessed March 2, 2011.
33. dcm4che2 DICOM Toolkit. Available at: <http://www.dcm4che.org/confluence/display/d2/dcm4che2+DICOM+Toolkit>. Accessed March 1, 2011.
34. Toolkit DCMTK-DICOM. Available at: <http://dicom.offis.de/dcm4che2.php.en>. Accessed March 1, 2011.
35. Highcharts JS. Available at: <http://highcharts.com/>. Accessed March 2, 2011.
36. Cook TS, Zimmerman SL, Steingall S, Boonn WW, Kim W. An algorithm for intelligent sorting of CT-related dose parameters. *Proc SPIE* 2011;7961:79612H.

Informatics in Radiology

RADIANCE: An Automated, Enterprise-wide Solution for Archiving and Reporting CT Radiation Dose Estimates

Tessa S. Cook, MD, PhD • Stefan L. Zimmerman, MD • Scott R. Steingall, RT • Andrew D. A. Maidment, PhD • Woojin Kim, MD • William W. Boonn, MD

RadioGraphics 2011; 31:1833–1846 • Published online 10.1148/rg.317115048 • Content Codes: CT IN QA

Page 1835

It is important to remember that $CTDI_{vol}$ is not a measure of patient dose, but rather represents the energy output of the scanner (24).

Page 1837

RADIANCE populates the structured report template with the data from its database and generates an RDSR that can be automatically transmitted to the ACR DIR.

The contents of scanner-generated RDSRs can also be imported into the database, enabling users to centralize their dose monitoring in cases in which there is a mixture of CT scanners, only some of which can produce structured reports.

Page 1837

RADIANCE is able to process dose sheets generated by the four major CT vendors—Siemens, GE, Toshiba, and Philips—as well as from the CereTom portable head CT scanner by NeuroLogica (Danvers, Mass).

Page 1840

RADIANCE has a small footprint and can run on any standard, up-to-date Windows personal computer (PC).

Page 1842

RADIANCE is compliant with the IHE REM profile as both a dose information consumer and a dose information reporter.